



CRT Winter 4x10ml Vial Blood Specimen Shipper

Model Number: ATS1-4V48H-15P

Description: Multi-use blood specimen shipper for 4x10ml vials with 12" x 8" x 8" external dimensions. Validated per ISTA 7D testing to keep payload above +15°C for 48 hours in Winter.



Holdover Time per ISTA 7D		
Test	Payload Temperature	Holdover Time
Summer Cycle	-	-
Winter Cycle	Above +15°C	48 hours

Product Specifications			
External Dimensions	12" x 8" x 8"	Payload Capacity	4 x 10mL vials
Actual Weight	5 lbs	Dimensional Weight	6 lbs
Case Qty	6	Pallet Qty	72
Estimated Number of Uses		12	

#	Bill of Materials	Weight	Qty	Conditioning Instructions
A	MRF16120808 PUF Shipper	1.2 lbs	1	Store at room temperature of +20°C to +25°C for at least 24 hours prior to use. Condition at room temperature of +20°C to +25°C for at least 24 hours prior to use. Store at room temperature of +20°C to +25°C for at least 24 hours prior to use.
B	OM18P-5X8-400 PCM CRT Pack	1.35 lbs	2	
C	RLC50-0/30 +30C Temperature Strip	NA	2	
D	S-8136Wx19 Velcro Self-Grip Strap	NA	1	
E	S-4129 Outer Box	0.5 lbs	1	
F	TMP (optional) Sensor	NA	1	